


Vedr.: Nyt program til M2/X2	PUP nr.:	20-011
	Udarbejdet af:	LPA
	Dato:	26-10-2020
	Intern <input type="checkbox"/>	Ekstern <input checked="" type="checkbox"/>

Det nye program til M2/X2 inkluderer en opdateret batteriindikator, så den viser afladningen mere lineært.

Programmet opdateres efter behov. Har man allerede mc2009_app_01.00.0033c er det ok, men ældre versioner vil være en god ide at opdatere.

mc2009_app_01.00.0033d

Find programmet på vores hjemmeside på produksiden under fanen Tekniker



Egenskaber Teknisk information Brochure/Manualer FAQ **Tekniker**

PROGRAMMERING

MC X1 15 KM PROGRAM	Download
MC X1 10 KM PROGRAM	Download
MC X2 15 KM PROGRAM	Download
MC X2 10 KM PROGRAM	Download
MC M2/X2 DISPLAY PROGRAM	Download
MC M2/X2 RESET SERVICE PIN	Download
PUP_20-010 M2/X2 Software rev 32	Download
PUP_20-008 Amputationsstøtte Ergo std	Download
PUP_20-006 M2/X2 brugerpanel	Download
PUP_20-005 Eberspächer timer	Download
PUP_20-004 Generel info X-model	Download
PUP_20-001 Spændingsudligner på MC	Download
mc2009_app_rev 33d	Download

Man skal have et tekniker-login for at kunne hente programmet.

Den kan også hentes på Intranettet under Produkter – Scootere

NYHETER | PERSONAL | PRODUKTER

Produkter > Scootere

mc2009_app_01.00.0033d



Klaus Smedegaard
R & D Manager




medemagruppen

Regarding: New program for M2/X2	PUP No.:	20-011
	Prepared by:	LPA
	Date:	26-10-2020
	Internal <input type="checkbox"/>	External <input checked="" type="checkbox"/>

The new program for M2/X2 include the updated battery indicator, so it shows the discharge in a more linear way.

Update program when needed. If your program is mc2009_app_01.00.0033c then it is ok, but older version would be a good idea to update.

mc2009_app_01.00.0033d

Find the program on our webside on the product page under the Technician tab.



Egenskaber Teknisk information Brochure/Manualer FAQ **Tekniker**

PROGRAMMERING	
MC X1 15 KM PROGRAM	Download
MC X1 10 KM PROGRAM	Download
MC X2 15 KM PROGRAM	Download
MC X2 10 KM PROGRAM	Download
MC M2/X2 DISPLAY PROGRAM	Download
MC M2/X2 RESET SERVICE PIN	Download
PUP_20-010 M2/X2 Software rev 32	Download
PUP_20-008 Amputationsstøtte Ergo std	Download
PUP_20-006 M2/X2 brugerpanel	Download
PUP_20-005 Eberspächer timer	Download
PUP_20-004 Generel info X-model	Download
PUP_20-001 Spændingsudligner på MC	Download
mc2009_app rev 33d	Download

To retrieve the program, you must have a password for technicians.



Klaus Smedegaard
R & D Manager



medemagruppen